



copy of Termoblocco MyBlock Mini con raffreddamento da -0°C a 100°C a blocchi intercambiabili da 0.2 a 50ml

Termoblocco MyBlock Mini con raffreddamento da -25°C a 100°C a blocchi intercambiabili da 0.2 a 50ml. Con coperchio trasparente. Senza blocchi da acquistare separatamente



Product codes:

Reference: 81 912 0121

Product gallery:



Product description:

Con la sua struttura compatta, l'ampia gamma di temperature e le opzioni di più blocchi intercambiabili, myBlock Mini è lo strumento perfetto per le incubazioni di laboratorio. I semplici controlli tramite touchpad e un display digitale consentono la selezione della temperatura "set and walk away" con un'accuratezza eccellente.

Perfetto per l'uso da banco, l'adattatore opzionale da 12 V espande l'utilità di myBlock Mini all'utilizzo. Con i suoi 11cm di larghezza, myBlock Mini è davvero il primo incubatore a blocchi intercambiabili personale.

Product features:



TBH - TYPE: heater
TBH - No BLOCKS: 1
TBH - WITH COOLING: no
TBH - TEMPERATURE MAXIMUM: +160°C
TBH - THE MINIMUM TEMPERATURE REACHED: 30° below ambient (10°C)
TBH - POWER HEATING: 300W
TBH - RESOLUTION SETTING: 0,1°C
TBH - PRECISION: 0,5°C
TBH - UNIFORMITY: 0,5°C
TBH - DISPLAY: LED
TBH - TIMER: da 1m a 99h59m
TBH - PROGRAMS: not provided
TBH - AUTONOMY: -
TBH - COMPLIANCE: CE
TBH - FREE: no
TBH - INDEPENDENT BLOCK: no
TBH - OPTIONAL ACCESSORIES: blocchi varie misure
TBH - OPTIONAL ACCESSORIES: blocks various sizes
TBH - POWER: 220 V 50/60 Hz
TBH - WARRANTY: 12 months on manufacturing defects
TBH - FEATURES: Termoblocco a singolo blocco intercambiabile con temperatura fino a 160°C
TBH - TO BE USED WITH: provette Ø40mm
TBH - TO BE USED WITH: provette Ø28mm
TBH - TO BE USED WITH: provette Ø26mm
TBH - TO BE USED WITH: provette Ø20mm
TBH - TO BE USED WITH: provette Ø19mm
TBH - TO BE USED WITH: provette Ø15mm
TBH - TO BE USED WITH: provette Ø10mm
TBH - TO BE USED WITH: provette Ø30mm
TBH - TO BE USED WITH: microtubi 5.0ml
TBH - TO BE USED WITH: provette Ø13mm
TBH - TO BE USED WITH: provette Ø12mm
TBH - TO BE USED WITH: test tubes, Ø16mm
TBH - TO BE USED WITH: microtubes 2.0 ml
TBH - TO BE USED WITH: microtubes 1.5 ml
TBH - TO BE USED WITH: microtubes 0.5 ml
TBH - TO BE USED WITH: microtubes 0.2 ml