




tips to gel with filter neutral from 0.1 to 20ul, sterile with filter rack

tips to gel with filter - Free from DNA, Dnase and Rnase - Apirogeni EC neutral from 0.1 to 20ul,
sterile with filter rack-96 pieces-Cf.960pz



Product codes:

Reference: 64 885 0261

Product gallery:



Product description:

Ferrules made of transparent, with a section of the capillary is long and thin for gel loading allow extreme accuracy in sample loading in applications such as DNA sequencing and lelettroforesi of proteins. The capillary thin, it is particularly suitable for handling the gel in the wells is extremely narrow. Tips for gel available in bags of 1,000 pieces non-sterile and in rack with lid of 96 pieces sterile and free from DNA, Dnase and Rnase Apirogeni. Compatible with pipettes Gilson Pipetman® P20-P100-P200 Eppendorf® Biohit® Proline & mLine Socorex® Acura HTL® Discovery Nichiryo® Thermo Scientific Finnpiette and the like. Also available with filter (barrier spray) to avoid contamination from the air and the cross-contamination of the samples. Material: polypropylene. Autoclavable.

Product features:

PNT - SIZE: nd

THE PNT - TYPE: gel loading

THE PNT - TYPE: with filter

PNT - CAPACITY MAX: 20ul

PNT - RANGE: 0.1-20ul

PNT - FILTER: caricamento gel

PNT - FORM: capillary



PNT - TIP: capillary round
PNT - STERILITY: you
PNT - PURITY: DNA / DNase / RNase / Pyrogen FREE
PNT - AUTOCLAVABLE: no
PNT - SIZED: no
PNT - COMPLIANCE: CE
PNT - APPLICATIONS: loading gel
PNT - MATERIAL: PP
PNT - COLOR: neutral
PNT - FEATURES: tip for gel loading by 20ul rack, sterile
PNT - COMPATIBLE WITH: Biohit mLine
PNT - COMPATIBLE WITH: Biohit Proline
PNT - COMPATIBLE WITH: Eppendorf
PNT - COMPATIBLE WITH: Finnpipette
PNT - COMPATIBLE WITH: Socorex Acura
PNT - COMPATIBLE WITH: Nichiryo
PNT - COMPATIBLE WITH: Gilson Pipetman
PNT - COMPATIBLE WITH: HTL Discovery
PNT - PACKAGING: on the rack