



Microtubi qPCR 0.2ml con tappo piatto ottico - ultra trasp

Cf 1000pz - Microtubi qPCR 0.2ml con tappo piatto ottico - ultra trasp. Cf 1000pz

Product codes:

Reference: 12 879 6741



Product gallery:



Product description:

Microtubi da 0.2ml con tappo ottico piatto anti-impronta. Compatibili con molti termociclatori per PCR e PCR Real Time che accettano plastica per biologia molecolare a profilo regolare. Come ABI / LT, Bio-Rad, Eppendorf e altri.

- Profilo regolare
- Ultra trasparenti
- Pareti sottili
- Miscela in PP specifico per applicazioni in PCR
- Basso grado di evaporazione
- Bassa adesione di DNA e proteine sulle pareti del tubo
- Tappo ottico
- DNA, RNA, DNase, RNase e Pyrogen FREE
- Autoclavabili

Product features:

MTB - CAP: plato
MTB - TIPO DE: PCR
MTB - VOLUMEN: 0.2ml
MTB - CAPACIDAD: 0.2ml



MTB - CIERRE: la presión
MTB - CAP FORABILE: nd
MTB - INSERTA LA TAPA: no
MTB - ÁREA ÓPTICA: 12.6mm
MTB - FORMULARIO: ajustar
MTB - FONDO: cónico
MTB - AUTO DE PIE: no
MTB - MATERIAL: PP Low Binding
MTB - TAMAÑO: no
MTB - OPACIDAD: ultra-transparente
MTB - ÁREA DE ESCRITURA: no
MTB - ETIQUETA: ausente
MTB - RANGO DE TEMPERATURA: -30°C a +100°C
MTB - TEMPERATURA MAX: 100°C
MTB - TEMPERATURA MÍNIMA: -30°C
MTB - PUREZA: DNA / RNA / DNase / RNase / Pyrogen FREE
MTB - ESTERILIZABLE EN AUTOCLAVE: usted
MTB - ESTERILIDAD: no
MTB - TAMAÑO: h 20.9/21.8mm
MTB - CARACTERÍSTICAS: microtubi da 0.2ml regular profile con tappo piatto per PCR Real-Time ultra-transparenti - raccomandate quando il segnale è ottenuto dal fondo del microtubo
MTB - APROPIADO PARA: preparación de la muestra
MTB - APROPIADO PARA: almacenamiento
MTB - APROPIADO PARA: PCR
MTB - APROPIADO PARA: PCR Real Time
MTB - APROPIADO PARA: biología molecular
MTB - RESISTENTE: la congelación a -20°C
MTB - EMBALAJE: sobres a granel
MTB - PAREDES: fina de 0,3 mm
MTB - LA EVAPORACIÓN DE LA MUESTRA: el menor de 3%
MTB - UNIÓN: Menor del ADN 0.3%
MTB - UNIÓN: las proteínas de menos de 0,8%
MTB - INTENSIDAD DE SEÑAL: 100
MTB - INTERFERENCIA DE LA SEÑAL: 100
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2400-24 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9600-96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Deep Well
MTB - COMPATIBLE CON: Bio-Rad Laboratories-C1000-96-Well Fast
MTB - COMPATIBLE CON: Bio-Rad Laboratories-C1000-Dual 48/48 Fast
MTB - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Deep Well
MTB - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-96-Well Fast
MTB - COMPATIBLE CON: Bio-Rad Laboratories-C1000 Touch-Dual 48/48 Fast
MTB - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Connect Real-Time PCR Detection System-96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-CFX96 Real-Time PCR Detection System-96-Deep Well
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-GeneCycler-24 x 0.2ml



MTB - COMPATIBLE CON: Bio-Rad Laboratories-iCycler-48 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-iCycler-96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-MJ Mini-48-12 48 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-MyCycler-96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-PTC-100-96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Deep Well
MTB - COMPATIBLE CON: Bio-Rad Laboratories-S1000-96-Well Fast
MTB - COMPATIBLE CON: Bio-Rad Laboratories-S1000-Dual 48/48 Fast
MTB - COMPATIBLE CON: Bio-Rad Laboratories-T100-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler-96-77 96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler ep gradient Silver-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler gradient-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler gradient-96-77 96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler nexus-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler nexus eco-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler nexus gradient-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler nexus gradient eco-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler nexus GSX1-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler nexus GSX1e-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler nexus SX1-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler nexus SX1e-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler Personal-25-16 25 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler pro-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler pro S-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Thermomixer C-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Thermomixer comfort-96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-100-96 x 0.2ml AgV
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-150 MiniCycler-25 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-30-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-48-48 Alpha unit 48 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Alpha unit 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-96v Moto Alpha Unit 96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler-96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-T1 Thermocycler Combi-96-77 96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler-48 x 0.2ml



MTB - COMPATIBLE CON: Whatman Biometra-T3 Thermocycler Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler-48 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-T3000 Thermocycler Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TGradient-96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TPersonal-48 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TPersonal Combi-48-18 48 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional-96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional Basic-96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional Basic Gradient-96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional Standard Gradient-96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO-48 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO combi-48-30 30 x 0.5ml
MTB - COMPATIBLE CON: Whatman Biometra-TProfessional TRIO combi-48-30 48 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TRobot-96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-UNO II-96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-UNO-thermoblock-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-7300 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-7500 Real-Time PCR System-96 x 0.2ml Standard
MTB - COMPATIBLE CON: Applied Biosystems-7900HT Fast Real-Time PCR System-96 x 0.2ml Standard
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2700-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 2720-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 5700-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700-Dual 96-Well 96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700 Aluminum-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-GeneAmp PCR System 9700 Gold-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-ProFlex™ PCR System-32 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-ProFlex™ PCR System-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 3 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-QuantStudio™ 5 Real-Time PCR System-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-SimpliAmp™-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-Veriti™ 96-Well Thermal Cycler-96 x 0.2ml
MTB - COMPATIBLE CON: Applied Biosystems-Veriti™ Dx 96-Well Thermal Cycler-96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine-Chromo4 96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad-Chromo4 96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Dyad Disciple-Chromo4 96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-iQ Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-iQ5 Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-MiniOpticon-48-12 48 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-MyiQ Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE CON: Bio-Rad Laboratories-MyiQ2 Real-Time PCR-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler ep realplex-96 x 0.2ml
MTB - COMPATIBLE CON: Eppendorf-Mastercycler ep realplex Silver-96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-200 DNA Engine-Chromo4 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-220 DNA Engine Dyad-Chromo4 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-221 DNA Engine Disciple-Chromo4 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-225 DNA Engine Tetrad-Chromo4 96 x 0.2ml
MTB - COMPATIBLE CON: MJ Research Inc.-PTC-240 DNA Engine Tetrad 2-Chromo4 96 x 0.2ml
MTB - COMPATIBLE CON: Stratagene-Mx3000P™-96 x 0.2ml
MTB - COMPATIBLE CON: Stratagene-Mx3005P™-96 x 0.2ml
MTB - COMPATIBLE CON: Stratagene-Mx4000™-96 x 0.2ml
MTB - COMPATIBLE CON: Stratagene-Robocycler Gradient 96-96 x 0.2ml
MTB - COMPATIBLE CON: Whatman Biometra-TOptical-96 x 0.2ml



MTB - COMPATIBLE CON: Bio-Rad Laboratories iQ Real-Time PCR 96 x 0.2ml

MTB - COMPATIBLE CON: Bio-Rad Laboratories iQ5 Real-Time PCR 96 x 0.2ml